



California State University of Bakersfield, Department of Chemistry

## Briggs-Rauscher Oscillating Color Change Reaction



### Standards:

5-PS1-4. Conduct an investigation to determine whether the mixing of two or more substances results in new substances.

### Introduction:

*The Briggs-Rauscher reaction, also known as 'the oscillating clock', is one of the most common demonstrations of a chemical oscillator reaction. The reaction begins when three colorless solutions are mixed together. The color of the resulting mixture will oscillate between clear, amber, and deep purple for about 3-5 minutes. The solution ends up as a purple-black mixture.*

### Materials:

- Distilled water
- 21.5 g potassium iodate
- 2.25 mL sulfuric acid
- 7.8 g malonic acid
- 1.7 g manganese sulfate monohydrate
- 2 g of vitex starch
- 200 mL of 30% hydrogen peroxide
- 3 600 mL beakers
- 2 L beaker

- Magnetic stir bar
- Stir plate

**Safety:**

- Always wear eye protection and gloves when doing chemistry experiments
- Conduct this experiment in a well-ventilated area.

**Procedure:**

## 1. Solution A:

Add 21.5 g potassium iodate ( $\text{KIO}_3$ ) to ~400 mL distilled water. Stir in 2.5 mL sulfuric acid ( $\text{H}_2\text{SO}_4$ ). Continue stirring until the potassium iodate is dissolved. Dilute to 500 mL.

## 2. Solution B:

Add 7.8 g malonic acid ( $\text{HOOCCH}_2\text{COOH}$ ) and 1.7 g manganese sulfate monohydrate ( $\text{MnSO}_4 \cdot \text{H}_2\text{O}$ ) to ~400 mL distilled water. Add 2 g of vitex starch. Stir until dissolved. Dilute to 500 mL.

## 3. Solution C:

Dilute 200 mL of 30% hydrogen peroxide ( $\text{H}_2\text{O}_2$ ) to 500 mL.

4. Place the stirring bar into the large beaker.
5. Pour solutions A and B into the beaker.
6. Turn on the stirring plate. Adjust the speed to produce a large vortex.
7. Add solution C into the beaker. Enjoy!

**Data and Observations:**

What did you see? Anything you were not expecting?

**References:**

Helmenstine, Ph. D, Anne Marie. Briggs-Rauscher Oscillating Color Change Reaction.chemistry.about.

<http://chemistry.about.com/cs/demonstrations/a/aa050204a.htm> (accessed July 31, 2013)